

Belimumab Biosimilar (Anti-BAFF) Antibody

Product Information

Product SKU:

HDBS0044

Size:

100 µg

Isotype:

IgG1

Host:

Homo sapiens

Category:

Biosimilar Antibody

Uniprot:

Q9Y275

Reactivity:

Human

Synonyms:

TNFSF13B,BAFF,BLYS,CD257,DTL,TALL1,THANK,TNFSF20,ZTNF4,TALL-1

Common Name:

LymphoStat-B, hBlySmAb-1.1

Target:

BAFF

Formulation:

Powder

Buffer:

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions.



Guarantee:

12 months from date of dispatch

Storage and Stability:

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature.

Applications

ELISA

Recommended Dilutions:

ELISA 1:5000-10000

Research Use Only:

Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

