

## Bevacizumab Biosimilar (Anti-VEGFA) Antibody

### Product Information

**Product SKU:**

HDBS0059

**Size:**

100 µg

**Isotype:**

IgG1

**Host:**

Humanized

**Category:**

Biosimilar Antibody

**Uniprot:**

P15692

**Reactivity:**

Human

**Synonyms:**

rhuMAb-VEGF

**Common Name:**

MVCD1,VEGF,VPF

**Target:**

VEGFA

**Formulation:**

Powder

**Buffer:**

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions.



**Guarantee:**

12 months from date of dispatch

**Storage and Stability:**

Store at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$  for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at  $-80^{\circ}\text{C}$  (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature.

**Applications**

ELISA

**Recommended Dilutions:**

ELISA 1:5000-10000

**Research Use Only:**

Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

