

Shionogi Anti-CCR8 Biosimilar Antibody

Product Information

Product SKU:

HDBS0061

Size:

100 µg

Isotype:

IgG1

Host:

Humanized

Category:

Biosimilar Antibody

Uniprot:

P51685

Reactivity:

Human

Synonyms:

CC-CKR-8, CCR-8, CDw198, CKRL1, CMKBR8, CMKBRL2, CY6, GPRCY6, TER1

Common Name:

盐野义 (Shionogi Inc.) 10A11

Target:

CCR8

Formulation:

Powder

Buffer:

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions.



Guarantee:

12 months from date of dispatch

Storage and Stability:

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature.

Applications

ELISA, Flow Cyt

Recommended Dilutions:

ELISA 1:5000-10000; Flow Cyt 1:100

Research Use Only:

Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

