

Prezalumab Biosimilar (Anti-B7-H2) Antibody

Product Information

Product SKU:

HDBS0075

Size:

100 µg

Isotype:

IgG2

Host:

Homo sapiens

Category:

Biosimilar Antibody

Uniprot:

O75144

Reactivity:

Human

Synonyms:

ICOSLG, B7-H2, B7H2, B7RP-1, B7RP1, CD275, GL50, ICOS-L, ICOSL, LICOS, ICOS ligand

Common Name:

AMG557, AMG-557, MEDI-5872

Target:

B7-H2

Formulation:

Powder

Buffer:

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions.



Guarantee:

12 months from date of dispatch

Storage and Stability:

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature.

Applications

Flow Cyt

Recommended Dilutions:

Flow Cyt 1:100

Research Use Only:

Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

