

## Plozalizumab Biosimilar (Anti-CCR2) Antibody

### Product Information

**Product SKU:**

HDBS0080

**Size:**

100 µg

**Isotype:**

IgG1

**Host:**

Humanized

**Category:**

Biosimilar Antibody

**Uniprot:**

P41597

**Reactivity:**

Human

**Synonyms:**

CC-CKR-2, CCR2A, CCR2B, CD192, CKR2, CKR2A, CKR2B, CMKBR2, MCP-1-R

**Common Name:**

MLN-1202, MLN1202, anti-CCR2, hu1D9

**Target:**

CCR2

**Formulation:**

Powder

**Buffer:**

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions.



**Guarantee:**

12 months from date of dispatch

**Storage and Stability:**

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature.

**Applications**

Flow Cyt

**Recommended Dilutions:**

Flow Cyt 1:100

**Research Use Only:**

Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

