

Anti-CD19 Biosimilar (FMC63) Antibody

Product Information

Product SKU:

HDBS0084

Size:

100 µg

Isotype:

IgG2a

Host:

Mus musculus

Category:

Biosimilar Antibody

Uniprot:

P15391

Reactivity:

Mouse

Synonyms:

CD19,B4,CVID3,MGC12802

Common Name:

FMC63

Target:

CD19

Formulation:

Powder

Buffer:

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions.



Guarantee:

12 months from date of dispatch

Storage and Stability:

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature.

Applications

Flow Cyt

Recommended Dilutions:

Flow Cyt 1:100

Research Use Only:

Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

