

Monalizumab Biosimilar (Anti-NKG2A) Antibody

Product Information

Product SKU:

HDBS0085

Size:

100 µg

Isotype:

IgG4

Host:

Humanized

Category:

Biosimilar Antibody

Uniprot:

P26715

Reactivity:

Human

Synonyms:

CD159A, NKG2, NKG2A

Common Name:

IPH-2201, IPH2201, NN-8765, NNC-0141-0000-0100, NNC-141-01000, anti-NKG2A, humZ270

Target:

NKG2A

Formulation:

Powder

Buffer:

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions.



Guarantee:

12 months from date of dispatch

Storage and Stability:

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature.

Applications

Flow Cyt

Recommended Dilutions:

Flow Cyt 1:100

Research Use Only:

Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

