

## Ulocuplumab Biosimilar (Anti-CXCR4) Antibody

### Product Information

**Product SKU:**

HDBS0089

**Size:**

100 µg

**Isotype:**

IgG4

**Host:**

Homo sapiens

**Category:**

Biosimilar Antibody

**Uniprot:**

P61073

**Reactivity:**

Human

**Synonyms:**

CD184, D2S201E, FB22, HM89, HSY3RR, LAP-3, LAP3, LCR1, LESTR, NPY3R, NPYR, NPYRL, NPYY3R, WHIM, WHIMS

**Common Name:**

BMS-936564, MDX-1338

**Target:**

CXCR4

**Formulation:**

Powder

**Buffer:**

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions.

**Guarantee:**

12 months from date of dispatch

**Storage and Stability:**

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature.

**Applications**

Flow Cyt

**Recommended Dilutions:**

Flow Cyt 1:100

**Research Use Only:**

Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.