

Anti-CD27 Biotinylated Recombinant Rabbit Monoclonal Antibody

Product Information

Product SKU:

HDLA022

Size:

100 µg

Isotype:

Rabbit IgG

Host:

Rabbit

Clonality:

Monoclonal

Clone:

DM59

Category:

Recombinant Antibody

Reactivity:

Human

Synonyms:

CD27, TNFRSF7, S152, T14, Tp55

Formulation:

Powder

Purification Method:

Purified from cell culture supernatant by affinity chromatography

Buffer:

1XPBS

Storage and Stability:

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature.

Applications

ELISA, Flow Cyt

Recommended Dilutions:

ELISA 1:5000-10000; Flow Cyt 1:100

Background:

The protein encoded by this gene is a member of the TNF-receptor superfamily. This receptor is required for generation and long-term maintenance of T cell immunity. It binds to ligand CD70, and plays a key role in regulating B-cell activation and immunoglobulin synthesis. This receptor transduces signals that lead to the activation of NF-kappaB and MAPK8/JNK. Adaptor proteins TRAF2 and TRAF5 have been shown to mediate the signaling process of this receptor. CD27-binding protein (SIVA), a proapoptotic protein, can bind to this receptor and is thought to play an important role in the apoptosis induced by this receptor.

Research Use Only:

Research Use Only