

# Anti-B7-2 Biotinylated Recombinant Rabbit Monoclonal Antibody

## Product Information

**Product SKU:**

HDLA046

**Size:**

100 µg

**Isotype:**

Rabbit IgG

**Host:**

Rabbit

**Clonality:**

Monoclonal

**Clone:**

DM84

**Category:**

Recombinant Antibody

**Reactivity:**

Human

**Synonyms:**

CD86, B7-2, B70, CD28LG2, LAB72, MGC34413

**Formulation:**

Powder

**Purification Method:**

Purified from cell culture supernatant by affinity chromatography

**Buffer:**

1XPBS

**Storage and Stability:**

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature.

**Applications**

ELISA, Flow Cyt

**Recommended Dilutions:**

ELISA 1:5000-10000; Flow Cyt 1:100

**Background:**

This gene encodes a type I membrane protein that is a member of the immunoglobulin superfamily. This protein is expressed by antigen-presenting cells, and it is the ligand for two proteins at the cell surface of T cells, CD28 antigen and cytotoxic T-lymphocyte-associated protein 4. Binding of this protein with CD28 antigen is a costimulatory signal for activation of the T-cell. Binding of this protein with cytotoxic T-lymphocyte-associated protein 4 negatively regulates T-cell activation and diminishes the immune response. Alternative splicing results in several transcript variants encoding different isoforms.

**Research Use Only:**

Research Use Only