

Anti-BTN3A1 Biotinylated Recombinant Rabbit Monoclonal Antibody

Product Information

Product SKU:

HDLA054

Size:

100 µg

Isotype:

Rabbit IgG

Host:

Rabbit

Clonality:

Monoclonal

Clone:

DM92

Category:

Recombinant Antibody

Reactivity:

Human

Synonyms:

BTN3A1, BTF5, CD277, BTN3.1, BT3.1

Formulation:

Powder

Purification Method:

Purified from cell culture supernatant by affinity chromatography

Buffer:

1XPBS

Storage and Stability:

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature.

Applications

ELISA, Flow Cyt

Recommended Dilutions:

ELISA 1:5000-10000; Flow Cyt 1:100

Background:

The butyrophilin (BTN) genes are a group of major histocompatibility complex (MHC)-associated genes that encode type I membrane proteins with 2 extracellular immunoglobulin (Ig) domains and an intracellular B30.2 (PRYSPRY) domain. Three subfamilies of human BTN genes are located in the MHC class I region: the single-copy BTN1A1 gene (MIM 601610) and the BTN2 (e.g., BTN2A1; MIM 613590) and BTN3 (e.g., BNT3A1) genes, which have undergone tandem duplication, resulting in 3 copies of each.

Research Use Only:

Research Use Only