

Anti-Human IL-4 In Vivo Antibody - Low Endotoxin

Product Information

Product SKU:

IVMB0003

Size:

1mg, 5mg, 25mg, 50mg, 100mg

Concentration:

≥ 2.0 mg/ml

Isotype:

Rat IgG1 κ

Host:

Rat

Clone:

MP4-25D2

Category:

Monoclonal Antibody

Reactivity:

Human

Synonyms:

BCGF, BCDF, B Cell Stimulating Factor, BSF-1

Specificity:

Clone MP4-25D2 recognizes human IL-4.

Formulation:

This monoclonal antibody is aseptically packaged and formulated in 0.01 M phosphate buffered saline (150 mM NaCl) PBS pH 7.2 - 7.4 with no carrier protein, potassium, calcium or preservatives added.

Endotoxin Level:

< 1.0 EU/mg as determined by the LAL method

Purity:

≥95% monomer by analytical SEC · >95% by SDS Page

Immunogen:

CHO-expressed, recombinant human IL-4

Storage and Stability:

Functional grade preclinical antibodies may be stored sterile as received at 2-8°C for up to one month. For longer term storage, aseptically aliquot in working volumes without diluting and store at -80°C. Avoid Repeated Freeze Thaw Cycles.

Product Preparation:

Functional grade preclinical antibodies are manufactured in an animal free facility using only In vitro protein free cell culture techniques and are purified by a multi-step process including the use of protein A or G to assure extremely low levels of endotoxins, leachable protein A or aggregates.

Applications

CyTOF®, IHC, N, WB

Applications & Recommended Usage:

ELISA: Each lot of this antibody is quality control tested by ELISA.