

Anti-Mouse IFNγ [XMG1.2] In Vivo Antibody - Low Endotoxin

Product Information

Product SKU: IVMB0022	
Size: 1mg, 5mg, 25mg, 50mg, 100mg	
Concentration: ≥ 5.0 mg/ml	
Isotype: Rat IgG1 κ	
Host: Rat SSS VGENIE	
Clone: XMG1.2	
Category: Monoclonal Antibody	2

Reactivity:

Mouse

Synonyms:

Immune Interferon, Type II Interferon, T Cell Interferon, MAF, IFNG, IFG, IFI

Specificity:

RatAnti-Mouse Interferon Gamma (IFN-γ) (Clone XMG1.2) recognizes an epitope on Mouse IFN-γ. This monoclonal antibody was purified using multi-step affinity chromatography methods such as Protein A or G depending on the species and isotype.

Formulation:



This monoclonal antibody is aseptically packaged and formulated in 0.01 M phosphate buffered saline (150 mM NaCl) PBS pH 7.2 - 7.4 with no carrier protein, potassium, calcium or preservatives added.

Endotoxin Level:

< 1.0 EU/mg as determined by the LAL method

Purity:

≥95% monomer by analytical SEC · >95% by SDS Page

Immunogen:

Purified Recombinant Mouse IFN-y (>98%)

Storage and Stability:

Functional grade preclinical antibodies may be stored sterile as received at 2-8°C for up to one month. For longer term storage, aseptically aliquot in working volumes without diluting and store at -80°C. Avoid Repeated Freeze Thaw Cycles.

Product Preparation:

Functional grade preclinical antibodies are manufactured in an animal free facility using only In vitro protein free cell culture techniques and are purified by a multi-step process including the use of protein A or G to assure extremely low levels of endotoxins, leachable protein A or aggregates.

Applications

CyTOF®, ELISA Cap, ELISPOT, ICFC, IHC FF, N, WB

Applications & Recommended Usage:

ELISA Cap XMG1.2 antibody is useful as a detection antibody for a sandwich ELISA with a concentration range of 0.5-2.0 µg/ml as recommended.