

Anti-Mouse IL-12 In Vivo Antibody - Low Endotoxin

Product Information

Product SKU:

IVMB0023

Size:

1mg, 5mg, 25mg, 50mg, 100mg

Concentration:

≥ 2.0 mg/ml

Isotype:

Rat IgG2a κ

Host:

Rat

Clone:

C17.8

Category:

Monoclonal Antibody

Reactivity:

Mouse

Synonyms:

NKSF, TSF, Maturation Factor, Cytotoxic Lymphocyte Maturation Factor (CLMF), P40

Specificity:

Mouse, IL-12/IL-23 p40 subunit (monomer, homodimer and heterodimer IL-12 p35/p40 or IL-23 p19/p40)

Formulation:

This monoclonal antibody is aseptically packaged and formulated in 0.01 M phosphate buffered saline (150 mM NaCl) PBS pH 7.2 - 7.4 with no carrier protein, potassium, calcium or preservatives added.

Endotoxin Level:

< 1.0 EU/mg as determined by the LAL method

Purity:

≥95% monomer by analytical SEC · >95% by SDS Page

Immunogen:

Purified CHO expressed Recombinant Mouse IL-12 p70

Storage and Stability:

Functional grade preclinical antibodies may be stored sterile as received at 2-8°C for up to one month. For longer term storage, aseptically aliquot in working volumes without diluting and store at -80°C. Avoid Repeated Freeze Thaw Cycles.

Product Preparation:

Functional grade preclinical antibodies are manufactured in an animal free facility using only In vitro protein free cell culture techniques and are purified by a multi-step process including the use of protein A or G to assure extremely low levels of endotoxins, leachable protein A or aggregates.

Applications

ELISA, ICFC, IP, N, WB

Applications & Recommended Usage: