

Anti-Mouse CD8a (Ly 2.2) In Vivo Antibody - Low Endotoxin

Product Information

Product SKU:

IVMB0070

Size:

1mg, 5mg, 25mg, 50mg, 100mg

Concentration:

≥ 5.0 mg/ml

Isotype:

Rat IgG2b

Host:

Rat

Clone:

2.43

Category:

Monoclonal Antibody

Reactivity:

Mouse

Synonyms:

Ly 2.2, Ly-2, Ly-35, Ly-B, Lyt-2, Lyt2

Specificity:

Clone 2.43 recognizes an epitope on mouse CD8a.

Formulation:

This monoclonal antibody is aseptically packaged and formulated in 0.01 M phosphate buffered saline (150 mM NaCl) PBS pH 7.2 - 7.4 with no carrier protein, potassium, calcium or preservatives added.

Endotoxin Level:

< 1.0 EU/mg as determined by the LAL method

Purity:

≥95% monomer by analytical SEC · >95% by SDS Page

Immunogen:

Mouse CTL clone L3

Storage and Stability:

Functional grade preclinical antibodies may be stored sterile as received at 2-8°C for up to one month. For longer term storage, aseptically aliquot in working volumes without diluting and store at -80°C. Avoid Repeated Freeze Thaw Cycles.

Product Preparation:

Functional grade preclinical antibodies are manufactured in an animal free facility using only In vitro protein free cell culture techniques and are purified by a multi-step process including the use of protein A or G to assure extremely low levels of endotoxins, leachable protein A or aggregates.

Applications

Depletion, FA, FC, ICC, IF Staining, IHC FFPE, IP

Applications & Recommended Usage:

FC The suggested concentration for this 2.43 antibody for staining cells in flow cytometry is ≤ 1.0 µg per 10⁶ cells in a volume of 100 µl. Titration of the reagent is recommended for optimal performance for each application.