

Anti-Mouse IL-1 β In Vivo Antibody - Low Endotoxin

Product Information

Product SKU:

IVMB0085

Size:

1mg, 5mg, 25mg, 50mg, 100mg

Concentration:

≥ 5.0 mg/ml

Isotype:

Armenian Hamster IgG

Host:

Armenian Hamster

Clone:

B122

Category:

Monoclonal Antibody

Reactivity:

Mouse

Synonyms:

IL1B, IL-1, IL1-Beta, IL1F2

Specificity:

Armenian Hamster Anti-Mouse IL-1 β (Clone B122) recognizes Mouse IL-1 β . This monoclonal antibody was purified using multi-step affinity chromatography methods such as Protein A or G depending on the species and isotype.

Formulation:

This monoclonal antibody is aseptically packaged and formulated in 0.01 M phosphate buffered saline (150 mM NaCl) PBS pH 7.2 - 7.4 with no carrier protein, potassium, calcium or preservatives added.

Endotoxin Level:

< 1.0 EU/mg as determined by the LAL method

Purity:

≥95% monomer by analytical SEC · >95% by SDS Page

Immunogen:

E. coli -expressed, recombinant mouse IL-1 β

Storage and Stability:

Functional grade preclinical antibodies may be stored sterile as received at 2-8°C for up to one month. For longer term storage, aseptically aliquot in working volumes without diluting and store at -80°C. Avoid Repeated Freeze Thaw Cycles.

Product Preparation:

Functional grade preclinical antibodies are manufactured in an animal free facility using only In vitro protein free cell culture techniques and are purified by a multi-step process including the use of protein A or G to assure extremely low levels of endotoxins, leachable protein A or aggregates.

Applications

IP, N, WB

Applications & Recommended Usage: