

Anti-Mouse Ly-6G/Ly-6C (Gr-1) In Vivo Antibody - Ultra Low Endotoxin

Product Information

Product SKU:

IVMB0092

Size:

1mg, 5mg, 25mg, 50mg, 100mg

Concentration:

≥ 5.0 mg/ml

Isotype:

Rat IgG2b

Host:

Rat

Clone:

RB6-8C5

Category:

Monoclonal Antibody

Reactivity:

Mouse

Synonyms:

Ly-6G/LY-6C, Gr-1, Myeloid

Specificity:

RatAnti-Mouse Granulocytes (Clone RB6-8C5) recognizes Mouse Granulocytes . This monoclonal antibody was purified using multi-step affinity chromatography methods such as Protein A or G depending on the species and isotype. The antibody was also test by PCR and IMPACT1 certified (For more information contact your Sales Rep).

Formulation:



This monoclonal antibody is aseptically packaged and formulated in 0.01 M phosphate buffered saline (150 mM NaCl) PBS pH 7.2 - 7.4 with no carrier protein, potassium, calcium or preservatives added.

Endotoxin Level:

<0.5 EU/mg as determined by the LAL method

Purity:

≥98% monomer by analytical SEC · >95% by SDS Page

Immunogen:

Raised against granulocytes of mouse origin

Storage and Stability:

Functional grade preclinical antibodies may be stored sterile as received at 2-8°C for up to one month. For longer term storage, aseptically aliquot in working volumes without diluting and store at -80°C. Avoid Repeated Freeze Thaw Cycles.

Product Preparation:

Functional grade preclinical antibodies are manufactured in an animal free facility using only In vitro protein free cell culture techniques and are purified by a multi-step process including the use of protein A or G to assure extremely low levels of endotoxins, leachable protein A or aggregates.

Applications

CyTOF®, Depletion, FC, IHC FF, IHC FFPE, IP, WB

Applications & Recommended Usage:

FC The suggested concentration for this RB6-8C5 antibody for staining cells in flow cytometry is ≤ 0.06 µg per 10⁶ cells in a volume of 100 µl or 100 µl of whole blood. Titration of the reagent is recommended for optimal performance for each application.