

Rat IgG2a Isotype Control - Low Endotoxin

Product Information

Product SKU:

IVMB0199

Size:

1mg, 5mg, 25mg, 50mg, 100mg

Concentration:

≥ 5.0 mg/ml

Isotype:

Rat IgG2a

Host:

Rat

Clone:

Sun Jan 01 2023 00:00:00 GMT+0000 (Greenwich Mean Time)

Category:

Isotype Control

Reactivity:

Synonyms:

Specificity:

This Rat IgG2a isotype control (anti-Trinitrophenol + KLH) antibody has been tested against selected species' cells and tissues to assure minimal cross reactivity.

Formulation:

This monoclonal antibody is aseptically packaged and formulated in 0.01 M phosphate buffered saline (150 mM NaCl) PBS pH 7.2 - 7.4 with no carrier protein, potassium, calcium or preservatives added.

Endotoxin Level:

< 1.0 EU/mg as determined by the LAL method

Purity:

≥95% monomer by analytical SEC · >95% by SDS Page

Immunogen:

Storage and Stability:

Functional grade preclinical antibodies may be stored sterile as received at 2-8°C for up to one month. For longer term storage, aseptically aliquot in working volumes without diluting and store at -80°C. Avoid Repeated Freeze Thaw Cycles.

Product Preparation:

Functional grade preclinical antibodies are manufactured in an animal free facility using only In vitro protein free cell culture techniques and are purified by a multi-step process including the use of protein A or G to assure extremely low levels of endotoxins, leachable protein A or aggregates.

Applications

Applications & Recommended Usage: