

SARS-CoV-2 Spike NTD [2146] Low Endotoxin Antibody (HRP)

Product Information

Product SKU:

IVMB0227

Size:

50ug

Concentration:

0.5 mg/ml

Isotype:

Human IgG1

Host:

Clone:

2146

Category:

BioSimilar Antibody

Reactivity:

SARS-CoV-2 · Virus

Synonyms:

COV2-2146, SARS-CoV2 Spike NTD, COVID-19, 2019-nCoV, Severe acute respiratory syndrome coronavirus 2, SARS-CoV2

Specificity:

Anti-SARS-CoV-2 Spike NTD-HRP, clone 2146, specifically targets an epitope on the SARS-CoV-2 spike protein N-terminal domain.

Formulation:

Endotoxin Level:

Purity:

Immunogen:

Storage and Stability:

This horseradish peroxidase conjugated monoclonal antibody is stable when stored at 2-8°C. Do not freeze.

Product Preparation:

Applications

ELISA, IHC

Applications & Recommended Usage:

ELISA, IHC

