

IVMB0244

Product Information

Product SKU:

IVMB0244

Size:

50ug

Category:

Polyclonal Antibody

Reactivity:

Human

Antibody Information

Synonyms:

HRG1-β1, HRG, Breast Cancer Cell Differentiation Factor p45, Acetylcholine Receptor Inducing Activity [ARIA], Glial Factor, Neu Differentiation Factor

Specificity:

Goat Anti-Human Neuregulin 1-β1 (NRG1-β1) recognizes Human NRG1-β1. This antigen affinity purified polyclonal antibody was purified using a proprietary chromatographic technique that includes covalently immobilizing the antigen proteins or peptides to agarose based beads. This purification method enhances , reduces nonspecific binding of extraneous IgG and provides you with the most reliable reagent available for your early discovery research.

Immunogen:

Purified Recombinant Human NRG1-β1 (Accession # AAA58639)

Storage and Stability:

The lyophilized, biotinylated antigen affinity purified polyclonal antibody can be stored desiccated at -20°C to -70°C for up to twelve months from date of receipt. The reconstituted bitoin conjugate can be stored for at least four weeks at 2-8°C. For long-term storage of the reconstituted conjugate, aseptically aliquot into working volumes and store at -20°C to -70°C in a manual defrost freezer. Avoid Repeated Freeze Thaw Cycles. No detectable loss of activity was observed after six months.