

Mouse Anti-PD-1 (CD279) Fc Reduced (LALA-PG) Recombinant Monoclonal Antibody

Product Information

Product SKU:

IVMB0280

Size:

1mg, 5mg, 25mg, 50mg, 100mg

Concentration:

≤1.0 EU/mg

Isotype:

Mouse IgG2a k

Host:

Mouse-Rat Chimeric

Clone:

RMP1-14-mAbMods (LALA-PG)

Category:

Recombinant Monoclonal Antibody

Reactivity:

Mouse

Synonyms:

Programmed Death-1, CD279, PD1, PDCD1

Specificity:

Anti-PD1 antibody (clone RMP1-14) recognizes an epitope on mouse PD-1.

Formulation:

This recombinant monoclonal antibody is aseptically packaged and formulated in 0.01 M phosphate buffered saline (150 mM NaCl) PBS pH 7.2 - 7.4 with no carrier protein, potassium, calcium or preservatives added. Due to inherent biochemical properties of antibodies, certain



products may be prone to precipitation over time. Precipitation may be removed by aseptic centrifugation and/or filtration.

Endotoxin Level:

<1.0 EU/μg

Purity:

>95% monomer by analytical SEC · >95% by SDS Page

Immunogen:

Mouse PD-1 transfected BHK cells

Storage and Stability:

Functional grade preclinical antibodies may be stored sterile as received at 2-8°C for up to one month. For longer term storage, aseptically aliquot in working volumes without diluting and store at -80°C. Avoid Repeated Freeze Thaw Cycles.

Product Preparation:

Functional grade preclinical antibodies are manufactured in an animal free facility using in vitro cell culture techniques and are purified by a multi-step process including the use of protein A or G to assure extremely low levels of endotoxins, leachable protein A or aggregates.

Applications

B, FA, In Vivo, WB

Applications & Recommended Usage:

B, FA, In Vivo, WB