

ACTB Monoclonal Antibody



MACO0006

Product Information

Size:

50ug

Reactivity:

Human, Rat, Mouse, Monkey, Dog,
Chicken, Hamster, Rabbit, Pig, Sheep

Source:

Mouse

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:5000

Protein Background:**Gene ID:**

ACTB

Uniprot

P60709

Synonyms:

ACTB; Actin, cytoplasmic 1; Beta-actin


Immunogen:

Synthetic Peptide

Storage:

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Product Images



Western blot analysis of HeLa (1), Rat brain (2), Mouse brain (3), (A 2) diluted at 1:10000.