

HSPA1L Monoclonal Antibody




MACO0016


Product Information

Size: 50ug	Protein Background:
Reactivity: Human, Mouse, Rat	Gene ID: HSPA1L
Source: Mouse	Uniprot P34931/P08107
Isotype: IgG	Synonyms: Heat shock 70 kDa protein 1-like/Heat shock 70 kDa protein 1A/1B
Applications: ELISA, WB, IHC, IF	Immunogen: Synthetic Peptide
Recommended dilutions: WB:1:500-1:5000, IHC:1:50-1:500, IF:1:50-1:200	Storage: PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Product Images



Western blot analysis of 1) HeLa, 2) Mouse Brain, diluted at 1:2000.



IF analysis of HeLa with antibody (Left) and DAPI (Right) diluted at 1:100.