

Myosin Heavy Chain Monoclonal Antibody



MACO0032

Product Information

Size:

50ug

Protein Background:

Reactivity:

Human, Rat, Mouse, Fruit, Fly, Nematode

Gene ID:

Myosin Heavy Chain

Source:

Mouse

Uniprot

Isotype:

IgG

Synonyms:

Applications:

ELISA, IF

Immunogen:

Synthetic Peptide


Recommended dilutions:

IF:1:100-1:200

Storage:

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Product Images



Immunofluorescence Staining of nematode tissue with MYH mouse mAb (11 2) diluted at 1:100. (Provide by Tsinghua University) .