

POU2F1 Monoclonal Antibody




MACO0063

Product Information

Size: 50ug	Protein Background:
Reactivity: Human	Gene ID: POU2F1
Source: Mouse	Uniprot P14859/P09086
Isotype: IgG	Synonyms: POU domain class 2 transcription factor
Applications: ELISA, WB	Immunogen: Synthetic Peptide
Recommended dilutions: WB:1:500-1:5000	Storage: PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Product Images



Western blot analysis of 1) HeLa, 2) Jurkat, 3) HepG2, diluted at 1:2000