

Di-methyl-Histone H3(K27) Monoclonal Antibody



MACO0072

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Mouse

Isotype:

IgG

Applications:

ELISA, WB, IP

Recommended dilutions:

WB:1:500-1:5000, IP:1:200-1:2000

Protein Background:**Gene ID:**

Di-methyl-Histone H3(K27)

Uniprot

P68431/Q71DI3/P84243

Synonyms:

HIST1H3A; H3FA; HIST1H3B; H3FL; HIST1H3C; H3FC; HIST1H3D; H3FB; HIST1H3E; H3FD; HIST1H3F; H3FI; HIST1H3G;


Immunogen:

Synthetic Peptide


Storage:

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Product Images



Western blot analysis of HeLa, diluted at 1) 1:1000 2) 1:3000



1) Input: HeLa Cell Lysate 2) IP product: IP dilute 1:200