

GFP Monoclonal Antibody



MACO0093

Product Information

Size:

50ug

Protein Background:**Reactivity:****Gene ID:**

GFP

Source:

Mouse

Uniprot**Isotype:**

IgG

Synonyms:**Applications:**

ELISA, WB, IP

Immunogen:

Recombinant Protein


Recommended dilutions:

WB:1:500-1:5000, IP:1:200-1:2000


Storage:

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Product Images



Western blot analysis of GFP transfected Hela, diluted at 1) 1:5000 2) 1:10000



1) Input (control) 2) IP products, antibody dilution 1:200