

mCherry Monoclonal Antibody



MACO0094

Product Information

Size:

50ug

Protein Background:**Reactivity:****Gene ID:**

mCherry

Source:**Uniprot**

Mouse

Isotype:**Synonyms:**

IgG

Applications:**Immunogen:**

ELISA, WB

Recombinant Protein

Recommended dilutions:**Storage:**

WB:1:500-1:5000

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Product Images



Western blot analysis of mCherry recombinant protein, diluted at 1) 1:5000 2) 1:10000