

# ACTB Monoclonal Antibody



MACO0103

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Rat, Mouse, Monkey, Dog,  
Chicken, Hamster, Rabbit

**Source:**

Mouse

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

WB:1:500-1:5000, IHC:1:50-1:500

**Protein Background:**

**Gene ID:**

ACTB

**Uniprot**

P60709

**Synonyms:**

ACTB; Actin, cytoplasmic 1; Beta-actin

**Immunogen:**


Synthetic Peptide

**Storage:**


PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

## Product Images

---



Western blot analysis of 1) 293T, 2) Mouse Kidney tissue, 3) Hela, 4) Rat Heart tissue, 5) Rat Brain tissue, diluted at 1:5000.



IHC Staining of Human ovary tissue, diluted at 1:200