

Plant actin Monoclonal Antibody



MAC00107

Product Information

Size:

50ug

Reactivity:

Various Plants

Source:

Mouse

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:5000

Protein Background:

Gene ID:

Plant actin

Uniprot

Synonyms:


Immunogen:

Recombinant Protein

Storage:

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Product Images



Western blot analysis of Arabidopsis, diluted at 1) 1:5000 2) 1:10000