

Strep Tag II Monoclonal Antibody



MAC00113

Product Information

Size:

50ug

Reactivity:

Source:

Mouse

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:5000

Protein Background:

Gene ID:

Strep Tag II

Uniprot

Synonyms:


Immunogen:

Recombinant Protein

Storage:

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Product Images



Western blot analysis of recombinant protein, mAb diluted at 1:5000.