

FN1 Monoclonal Antibody



MACO0115

Product Information

Size:

50ug

Reactivity:

Mouse, Rat, Human

Source:

Mouse

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:5000

Protein Background:**Gene ID:**

FN1

Uniprot

P02751

Synonyms:

FN1; FN; Fibronectin; FN; Cold-insoluble globulin; CIG

Immunogen:

Synthetic Peptide

Storage:

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Product Images



Western blot analysis of 1) Mouse Brain Tissue, 2) Rat Brain tissue, diluted at 1:2000.