

# Histone H2B Monoclonal Antibody



MACO0119

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Mouse

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:500-1:5000

**Protein Background:**

**Gene ID:**

Histone H2B

**Uniprot**

P33778

**Synonyms:**

HIST1H2BB; H2BFF; Histone H2B type 1-B; Histone H2B.1; Histone H2B. f; H2B/f

**Immunogen:**


Synthetic Peptide

**Storage:**

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

## Product Images

---



Western blot analysis of 1) HeLa, 2) 3T3, 3) Raw264.7, 4) Rat Brain, 5) Rat Kidney, diluted at 1:2000.