

IgG (Fc specific, Heavy Chain) Monoclonal Antibody

MAC00124



Product Information

Size:

50ug

Reactivity:

Rabbit

Source:

Mouse

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:5000

Protein Background:**Gene ID:**

IgG (Fc specific, Heavy Chain)


Uniprot**Synonyms:****Immunogen:**

Recombinant Protein

Storage:

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Product Images



Western blot analysis of Reduced Rabbit IgG, H+L, diluted at 1:10000.