

# IgG (Light Chain Specific) Monoclonal Antibody



MACO0156

---

## Product Information

**Size:**

50ug

**Reactivity:**

Rabbit

**Source:**

Mouse

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:500-1:5000

**Protein Background:**

**Gene ID:**

IgG (Light Chain Specific)

**Uniprot**

**Synonyms:**

**Immunogen:**


Recombinant Protein

**Storage:**

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

## Product Images

---



Western blot analysis of Reduced Rabbit IgG, H+L, diluted at 1:10000.