

AQP4 Monoclonal Antibody



MACO0170

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Mouse

Isotype:

IgG1

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

WB:1:500-1:5000, IHC:1:50-1:500, IF:1:50-1:200

Protein Background:

Aquaporin 4 belongs to the Aquaporin family of integral membrane proteins that conduct water through the cell membrane.

Gene ID:

AQP4

Uniprot

P55087

Synonyms:

AQP-4; AQP4; aquaporin 4; aquaporin type4; Aquaporin-4; HMIWC2

Immunogen:

Recombinant Protein

Storage:

Mouse IgG1 in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images

N/A

N/A