

# ARRB1 Monoclonal Antibody



MAC00171

---

## Product Information

**Size:**

100ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Mouse

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:500-1:5000

**Protein Background:**

**Gene ID:**

ARRB1

**Uniprot**

P49407

**Synonyms:**

arrestin-beta-1

**Immunogen:**

Recombinant fragment, corresponding to amino acids of human ARRB1

**Storage:**

Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol

## Product Images

---

N/A

N/A