

# BECN1 Monoclonal Antibody



MACO0175

---

## Product Information

**Size:**

100ul

**Reactivity:**

Human, Rat, Mouse, Pig, Rabbit

**Source:**

Mouse

**Isotype:**

IgG1

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:500-1:5000

**Protein Background:**

Plays a central role in autophagy. Required for the abscission step in cytokinesis. May play a role in antiviral host defense. Protects against infection by a neurovirulent strain of Sindbis virus.

**Gene ID:**

BECN1

**Uniprot**

Q14457

**Synonyms:**

ATG6, &alpha;v&beta;30, beclin1

**Immunogen:**

Purified recombinant human BECN1 protein fragments expressed in E.coli.

**Storage:**

Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol

## Product Images

---

N/A

N/A