

# DsbA Monoclonal Antibody



MAC00218

---

## Product Information

**Size:**

100ul

**Reactivity:**

E.coli

**Source:**

Mouse

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:500-1:5000

**Protein Background:**

**Gene ID:**

DsbA

**Uniprot**

P0AEG4

**Synonyms:**

hypothetical protein

**Immunogen:**

Full Length (Escherichia coli) Disulfide Oxidoreductase, isolated from E. coli.

**Storage:**

Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol

## Product Images

---

N/A

N/A