

# IgG1 Monoclonal Antibody



MAC00319

---

## Product Information

**Size:**

100ul

**Reactivity:**

Human

**Source:**

Mouse

**Isotype:**

IgG1

**Applications:**

ELISA, IF, FC

**Recommended dilutions:**

IF:1:100-1:200

**Protein Background:**

This antibody can be used as a mouse IgG1 isotype control in flow cytometry and other applications.

**Gene ID:**

IgG1

**Uniprot**

**Synonyms:**

**Immunogen:**

Recombinant Protein

**Storage:**

Phosphate-buffered solution, pH 7.4, containing 0.09% sodium azide and 0.2% (w/v) BSA

## Product Images

---

N/A

N/A