

# NUP98 Monoclonal Antibody



MACO0328

---

## Product Information

**Size:**

100ul

**Reactivity:**

Human, Dog, Mouse

**Source:**

Mouse

**Isotype:**

IgG1

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:500-1:5000

**Protein Background:**

Nup98 and Nup96 play a role in the bidirectional transport across the nucleoporin complex (NPC). The repeat domain in Nup98 has a direct role in the transport.

**Gene ID:**

NUP98

**Uniprot**

P52948

**Synonyms:**

ADIR2, NUP96, NUP196

**Immunogen:**

Purified recombinant human NUP98 protein fragments expressed in E.coli.

**Storage:**

Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol

## Product Images

---

N/A

N/A