

TP53 Monoclonal Antibody



MACO0330

Product Information

Size:

100ul

Reactivity:

Human, Dog, Mouse, Rat

Source:

Mouse

Isotype:

IgG1

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:5000

Protein Background:

Gene ID:

TP53

Uniprot

P04637

Synonyms:

cellular tumor antigen p53 antigen NY-CO-13 phosphoprotein p53 p53 tumor suppressor transformation-related protein 53

Immunogen:

Recombinant fragment of human TP53 (C-terminus)

Storage:

Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol

Product Images

N/A

N/A