

Ubiquitin Monoclonal Antibody



MACO0394

Product Information

Size:

50ug

Protein Background:**Reactivity:**

Human, Mouse, Rat

Gene ID:

Ubiquitin

Source:

Mouse

Uniprot**Isotype:**

IgG

Synonyms:**Applications:**

ELISA, WB, IHC

Immunogen:

Synthetic Peptide of Ubiquitin


Recommended dilutions:

WB:1:500-1:5000, IHC:1:50-1:500


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.


Product Images



Immunohistochemical analysis of paraffin-embedded Human Stomach Carcinoma Tissue using Ubiquitin Mouse mAb diluted at 1:200.



Western blot analysis of 1) Hela Cell Lysate, 2) 3T3 Cell Lysate, 3) Rat Brain Tissue Lysate using Ubiquitin Mouse mAb diluted at 1:1000.



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma Tissue using Ubiquitin Mouse mAb diluted at 1:200.