

Luciferase Monoclonal Antibody



MAC00404

Product Information

Size:

50ug

Reactivity:

Firefly

Source:

Mouse

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:5000

Protein Background:

Gene ID:

Luciferase

Uniprot

P08659

Synonyms:


Immunogen:

Synthetic Peptide of Luciferase

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western blot analysis of Firefly Luciferase Transfected Hela Cell Lysate using Luciferase Mouse mAb diluted at 1:2000.