

TTR Monoclonal Antibody



MAC00442

Product Information

Size:

50ug

Reactivity:

Human

Source:

Mouse

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:5000, IHC:1:50-1:500

Protein Background:

Gene ID:

TTR

Uniprot

P02766

Synonyms:

Transthyretin (ATTR) (Prealbumin) (TBPA)


Immunogen:

Recombinant Protein


Storage:

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.


Product Images



Immunohistochemical analysis of paraffin-embedded Human Liver Carcinoma Tissue using TTR Mouse mAb diluted at 1:2000



Western blot analysis of Human Serum using TTR Mouse mAb diluted at 1:2000



Immunohistochemical analysis of paraffin-embedded Human Liver Carcinoma Tissue using TTR Mouse mAb diluted at 1:200.