

CALCA Monoclonal Antibody



MACO0553

Product Information

Size:

50ug

Reactivity:

Human

Source:

Mouse

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:50-1:500

Protein Background:

Gene ID:

CALCA

Uniprot

P01258

Synonyms:

PCT


Immunogen:

Recombinant Human PCT

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemical of paraffin-embedded human thyroid tissues using MACO0553 at dilution of 1:200