

Corticosterone Monoclonal Antibody



MACO0582

Product Information

Size:

50ug

Reactivity:

Source:

Mouse

Isotype:

IgG1

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

Corticosterone

Uniprot

Synonyms:

11-dehydro-17-hydroxycorticosterone, cortisone

Immunogen:

Corticosterone-BSA conjugate

Storage:

Preservative: 0.1mol/L NaCl, 15mmol/L NaN₃, PH7.2

Product Images

N/A

N/A