

Cholic acid Monoclonal Antibody



MACO0586

Product Information

Size:

50ug

Reactivity:

Source:

Mouse

Isotype:

IgG1

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

Cholic acid

Uniprot

Synonyms:

Cholate, cholanic acid, cholalin, colalin, cholalic acid, 3alpha,7alpha,12alpha-trihydroxy-5beta-cholanicacid, cholsaeure

Immunogen:

Cholic acid-BSA conjugate

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A