

The metabolite of furazolidone Monoclonal Antibody

MACO0587



Product Information

Size:

50ug

Reactivity:**Source:**

Mouse

Isotype:

IgG1

Applications:

ELISA

Recommended dilutions:**Protein Background:****Gene ID:**

The metabolite of furazolidone

Uniprot**Synonyms:**

CPAOZ

Immunogen:

The metabolite of furazolidone-BSA conjugate

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A