

Chenodeoxycholic acid Monoclonal Antibody



MACO0589

Product Information

Size:

50ug

Reactivity:

Source:

Mouse

Isotype:

IgG2b

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

Chenodeoxycholic acid

Uniprot

Synonyms:

Chenodiol, chenix, chenic acid, chenodeoxycholate, gallodesoxycholic acid, chendol, chenofalk, chenodesoxycholic acid, CDCA

Immunogen:

Chenodeoxycholic acid-BSA conjugate

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A